



Reference: PTXBC030011

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HsT19680; synaptic vesicle glycoprotein 2B; synaptic vesicle protein 2B homolog		
For research use only.			

GENE DESCRIPTION	synaptic vesicle glycoprotein 2B
GENE ID	9899
ORIGIN SPECIES	Human
ACCESSIONS	BC030011
VECTOR	pDONR223
NCBI SYMBOL	SV2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]