

HORMAD2-HORMA domain containing 2 Gene

Reference: PTXBC030013

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT46.2; HORMA domain-containing protein 2; HORMA domain containing 2		
For research use only.			

GENE DESCRIPTION	HORMA domain containing 2
GENE ID	150280
ORIGIN SPECIES	Human
ACCESSIONS	BC030013
VECTOR	pDONR223
NCBI SYMBOL	HORMAD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccactgctcagctttctcactgcatcacatacacaaggtcttaaggaaacagtttccatcccaatcactaatgagcatgagtcattgaaatgggaagaaacttttgcctactctcatgtataacatacctaagggcctgtttccagagagctcttatggagaacgccattg
gatgacctcagtttaaaatcctccgagaagataaaaaatgtccgggtcactgcatattatcagatggattcaaggtgtttgatgctttgaaaagagatacctacgtatggcagtgactgacatttacacagatcccatgggatctgagaagggtgactgagatgtaccagttcaaatcaaat
acacgaaagaaggagccactatggatttgacagtcatagcagcagtcacaagctttgaaagtggacaacaacatgaagataaagaaagccagtgcttctactgacccgtaaatgtatatactgacagagccttgagccacttctataatgtgtacttactgaaactccactactataat
gcagtgaccccatgattaccaaccctcgggttttaagaagggttaattcacactcctgctgtttgacaaggagcctatcaacgtgcaggtgggatttgcctcactggcttcatagcatgaaagttaaagtcagacagagggtacaaaagtgtattgtttggagaacaactctgttcgg
gagaacagcactactgagatcgccatcagggctctagactgtgatgaggaagaagaatgcaatgaccatattcaagaatgatttgtgtgcagtcagcaagttctgagtgctccaggaagaagaggaaggtcagtgaccagtgaaaggtcttccctaacagaaaatga