

TMEM196-transmembrane protein 196 Gene

Reference: PTXBC030104

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 196		
For research use only.			

GENE DESCRIPTION	transmembrane protein 196
GENE ID	256130
ORIGIN SPECIES	Human
ACCESSIONS	BC030104
VECTOR	pDONR223
NCBI SYMBOL	TMEM196
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcctctttcagcctgctgtatctgtggactattggggcatcctgaattttcagttctccgggcagtcacaaagaaaacttctccctataccactgcacctgcctccatgtctctcgctgcattgggatcggggcgtgcaactctctctctggtcacttgtcgactagccagttatg
aacagaggaggatgtttcagaaaggagcattccctgcctcactctcatgaaatggctgagaaaaatgagggtattga