

MGC42157-hypothetical locus MGC42157 Gene

Reference: PTXBC030111

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical locus MGC42157
GENE ID	439933
ORIGIN SPECIES	Human
ACCESSIONS	BC030111
VECTOR	pDONR223
NCBI SYMBOL	MGC42157
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaaggctaaactttttctccacatgctcacttcctcgtttgtcgtagctgctgcctcgaagggttaaagtaagggaaggagagattgagaggaaggacagctgtgaattgggtctcaccgtgtttccaggctgaagtacagtggttgatcagttcattgcaact
tctgtctctgagctcaagcaatcctcttgcctcagcctcgaagtagctgggactgactacaggcatgaaccaccacacctggcctaattttttgtatgttgcgaaaagacagtttcgccatgtgtacagaactactctcaaatcactattcacatcctggtag