

FAM171A2-family with sequence similarity 171, member A2 Gene

Reference: PTXBC030200

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM171A2; family with sequence similarity 171 member A2		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 171, member A2
GENE ID	284069
ORIGIN SPECIES	Human
ACCESSIONS	BC030200
VECTOR	pDONR223
NCBI SYMBOL	FAM171A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtttgggaaccggactctgtggcagctggcaccacagactcagaggggtgtggccaccctgccctcagttatcgcttgggcacotgggtgctgtcactgctgccgcctggcttctcaccactctgtgccttggcgtgttgacaagctgcccttgcgtctgtcagcctta
cctgtcctctgagcggccggccacgctcatctctatgaggacctggtgcacattctcctag