

## CD79B-CD79b molecule, immunoglobulin-associated beta Gene

Reference: PTXBC030210

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CD79b molecule; CD79b molecule, immunoglobulin-associated beta; CD79b antigen (immunoglobulin-associated beta); AGM6; B29; IGB; B-cell antigen receptor complex-associated protein beta chain; B-cell-specific glycoprotein B29; Ig-beta; immunoglobulin-associ |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | CD79b molecule, immunoglobulin-associated beta                                                    |
| <b>GENE ID</b>          | 974                                                                                               |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC030210                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CD79B                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggccaggctggcgttctctctgtgccagccactggatggcgttctgtctgtctcaggtacagaacccacgacgagcctgtggggtttgtctctatctccagctgtctggccctgcggctctgctctctgtctccaccctctctctgtctctactctttaccgctgtctcctctcatggccctgggacctgtgtgggtgtcagctgcacgtggtga