

SLC22A18AS-solute carrier family 22 (organic cation transporter), member 18 antisense Gene

Reference: PTXBC030237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BWR1B; BWSCR1B; ORCTL2S; SLC22A1LS; p27-BWR1B; beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein; Beckwith-Wiedemann region 1B; Beckwith-Wiedemann syndrome chromosome region 1, candidate b; organic cation transporter-like 2 antisen		

For research use only.

GENE DESCRIPTION	solute carrier family 22 (organic cation transporter), member 18 antisense
GENE ID	5003
ORIGIN SPECIES	Human
ACCESSIONS	BC030237
VECTOR	pDONR223
NCBI SYMBOL	SLC22A18AS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgggtgagaggggtgctggtgttctctctgagccccgcagggctgcagctccatctggccagcagagggcgagaaagagctgctgggcagctcggacgtgaccgctggaagtgtgtacagcgttctgacaacgttccggccaaaacgggtccgcccac
ttgtgtcaagataactgaaatgtcttctgtgtcctccgaggagaatgaaaggctgtggatgtagtgccttctctactgctctcagctggggcagaatgacaggacgtgaacatgcacaaaggccgggaactacgacccccgacagcgttccccgccccca
ccccgcgttcagggatttaggcactagcacagcacagtcacccctgggactagcaccctccgcagttctgctctgcagcaacatag