

GAL-galanin prepropeptide Gene

Reference: PTXBC030241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAL-GMAP; ETL8; GALN; GLNN; GMAP; galanin peptides; galanin prepropeptide; galanin-message-associated peptide; galanin-related peptide; galanin/GMAP prepropeptide; galanin and GMAP prepropeptide		

For research use only.

GENE DESCRIPTION	galanin prepropeptide
GENE ID	51083
ORIGIN SPECIES	Human
ACCESSIONS	BC030241
VECTOR	pDONR223
NCBI SYMBOL	GAL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaggcagcgccctctctctcctcctcctcgcgcggcccttctgcctctcggggctctgtgcggccaaggaaaaacgaggctggaccctgaacagcgcgggtacgtgctggccacatgccgttggaaccacaggtcattcagcgacaagaatggcc
tcaccagcaagcgggagctgcggcccgaagatgacatgaaaccagggaagctttgacaggtccatcctgaaacaatatcatgcgcacaatcattgagtttctgtcttttgcattcctcaagagccgggtgccctcgaccgctctggtatctccccgccgagctcctcagaagacatc
gagcggctcctga