

FAM30A-family with sequence similarity 30, member A Gene

Reference: PTXBC030286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM30A; C14orf110		

For research use only.

GENE DESCRIPTION	family with sequence similarity 30, member A
GENE ID	29064
ORIGIN SPECIES	Human
ACCESSIONS	BC030286
VECTOR	pDONR223
NCBI SYMBOL	FAM30A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagacccctctccaggaggaggagcctggagacacaggccccatctccaatggggacacttcaggagcggctctcaggtcccgagaagaccttctgccacaggagacacacggacatcggaacggacggaggaggatgtgcagttgcagcctttcagcaga
cgccctgagaacggagggtcaagagttggagcagacgacagggcctcagacctaa