



Reference: PTXBC030507

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ETS-related transcription factor Elf-1; E74-like factor 1 (ets domain transcription factor); E74 like ETS transcription factor 1		

For research use only.

GENE DESCRIPTION	E74-like factor 1 (ets domain transcription factor)
GENE ID	1997
ORIGIN SPECIES	Human
ACCESSIONS	BC030507
VECTOR	pDONR223
NCBI SYMBOL	ELF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]