

## TDRD3-tudor domain containing 3 Gene

Reference: PTXBC030514

|                               |                                                              |                     |                |
|-------------------------------|--------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                          | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | tudor domain-containing protein 3; tudor domain containing 3 |                     |                |
| <b>For research use only.</b> |                                                              |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tudor domain containing 3                                                                        |
| <b>GENE ID</b>          | 81550                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC030514                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TDRD3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgctgcattacagatgactgatgtcatataagttgcacagcagtagaatttagttatgtcaaaaaataagcctgaacacaccactggaactaaagttaagctctcagggcattgttgacataaaaaatggattcctgctcttgatgacttaaccacacagttcttggtggaagtggaaac  
accttattgagaatgggagttacagagaagcttatcaaaacacaatagaagcaatattggaactgaaggtggaccaccgcttttggcttttgacagaagtggtatctcatgtccaagtgatagcagagaacttgatcgaagaaaaacattgcaagttacaatgctgtcaaacctaca  
aatgataatgatgaatttgaagcaaaagacggcgtctatttgcgaagttgcaaaagacgaaggaacccaagacatttggaggaggtggtggtggtgctagaagtaattcacaatgaatgctgctgtaaccgaaatagggaagtttacagaaaaaagtcaacaaatcagaggga  
aacatgaaggtgtctatagagaactgtgttgatgagaaagctctgaagcacataacggaatgggtcctagtaaggaagcatcgaggcaagctctatggataatggcaacaacttagaagcagcactgaacgtactcttacaagcaataaacagaacacgtgtatgggtcctcctgagag  
gtagaggaaaaagcagggggcgaataagatctgaagatgaagggacctgggaatgcaaggccatcagcaccacacacattattgatttctggaatctaaaatgggaactttgaatgtggaagaacctaaatcacagccacagcagctcatcagggaacatcacagatcaaaatc  
tgagcaaaatggagtaaaagataaatcatctgagacatcctcctcgaatgataccagcagcgaagaaatgaaaaacccctcgttttcaaaagagactcccaaaatcaaaagtcagttttagaaggcagtggtgattacctaagaatagaggttctgaaagaccaagtcttctcagatctg  
aagtatgggtggaagacagaatcaaatgtgatagaccgtattctagatagacagaactaaagatacttcatactcttttaggttctcagcatagtgatggtgcttttaaaaaaagagataaactctatgcaagcagatcaggaaggtgctcctcttgcagaggcgaagaaaaatccactcctc  
aaggatctgtagattataataatcaaaacgtggaaaaagagaagccaaacatcttctgattttttatgacaggaaatcacaaacaataataatgaagcttctcagtggtataaaaatgaaaaacatttgaatgtaatactgattatcagaatccagttcgaaatgaatgatttctggtt  
ccaaatggagaagtagaataatgccactgaaaggagacgaatagaccattatgaagccagcagcagctgtcacagctgtaccctgtgatgataaaatatttacaatagtgggcccaaacgaagatctggccaattaaagccagaaaaatactagaatcatctatctctatggagtgcacaaa  
atgtgaaacctggagatgaatgtttgcaactttatgggaagacaacaagttttaccgggcagaaagttgaagccctccattctcgggtatgacacagctgtgtaaatctatgactacggaaactatgaagaggtgctactgagcaatatcaagccattcaaacagagcatgggaggaag  
aaggcacctacgatcaactctggagttccgtaggggaggtgatggccagcaagacgatccactcggccaaccaacagttttaccaaccacccgggtcgcgaactaa