



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MRPS27-mitochondrial ribosomal protein S27 Gene

Reference: PTXBC030521

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP-S27; S27mt; 28S ribosomal protein S27, mitochondrial; mitochondrial 28S ribosomal protein S27; mitochondrial ribosomal protein S27		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S27
GENE ID	23107
ORIGIN SPECIES	Human
ACCESSIONS	BC030521
VECTOR	pDONR223
NCBI SYMBOL	MRPS27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctctgatatggaaccaggctactctgcagagacccctcaagtgatggagaagaatgctgctcctcccaagaagacataaagctgtgtagaagaagcgtctgatgtctgggtgcagtgctgaaggctctgacttcagctgatgggcttcagaggagcagctcccaaatgatgaagac
aaccaagcctgcagaataaactctggagacagcttagatcactgcaggaagaacagagactccaagctctcaataactggaaacgalltaaggccttcatcttgaagctccaagctctggggcaaaattatgcagaagctctttaaagtctccaccagcctgtcgaaggaataactccacctgt
ggcaagaagacactccaccattatgcacaatactccacagctctcactatgcactgtatgacttgcctcagaagaacacagcaagcagaagcgaacgaacaggatgaccacgtcgaagaagaacccaagaagctcactctccttaa