



Reference: PTXBC030530

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DER3; HRD1; E3 ubiquitin-protein ligase synoviolin; HMG-coA reductase degradation 1 homolog; synovial apoptosis inhibitor 1, synoviolin; synoviolin 1		

For research use only.

GENE DESCRIPTION	synovial apoptosis inhibitor 1, synoviolin
GENE ID	84447
ORIGIN SPECIES	Human
ACCESSIONS	BC030530
VECTOR	pDONR223
NCBI SYMBOL	SYVN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]