



Reference: PTXBC030531

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf24; GIG39; NIBAN; protein Niban; cell growth-inhibiting gene 39 protein; family with sequence similarity 129 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 129, member A
GENE ID	116496
ORIGIN SPECIES	Human
ACCESSIONS	BC030531
VECTOR	pDONR223
NCBI SYMBOL	FAM129A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]