

CD99L2-CD99 molecule-like 2 Gene

Reference: PTXBC030536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD99B; MIC2L1; CD99 antigen-like protein 2; MIC2 like 1; MIC2-like protein 1; CD99 molecule like 2		

For research use only.

GENE DESCRIPTION	CD99 molecule-like 2
GENE ID	83692
ORIGIN SPECIES	Human
ACCESSIONS	BC030536
VECTOR	pDONR223
NCBI SYMBOL	CD99L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggcctggcgtcggcgttcctgtctgcctcgtttctccttggccaccctggtccagcgaggatctggggacttggatgattttaacctggaggatgcagtgaagaaacttctcagtaagcagccatgggaccacaccaccaccacacaaccaataggccagggaaccaccag
agctccggcaaaacctccaggttagtggaattggacttggctgatcttggatgatcaagatgatggccgaggaaccgggtataggaggagagagatggaacctgtaaccaccagaccaagaggccagtaaccaccagagctccagcaatactttaggaaatgatttgactt
ggctgatccctggatgatcgaatgatcgagatgatggccgaggaaccaattgctggaggaggaggttttcagacaaggatcttgaagacatagtaggggtggagaatacaaacctgacaagggttaagggtgatggccgtgacggcagcaatgacgacctggatctggcatg
gtggcagagcctggcaccattgccgggtggccagcgccctggccatggccctcatcgggtgctcctcagctacatctctaccagcagaagaagtctctcagcattcagcaggggtctcaacgcagactacgtgaaggagagaacctggaagccgtggtatgtgaggaacccc
aagtgaatactccacgttgcacacgcagctgcagagccgccgccgccgaaccagcccgatctga