

ARL17-ADP-ribosylation factor-like 17 Gene

Reference: PTXBC030570

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARL17; ARL17A; ARL17P1; ADP-ribosylation factor-like 17 pseudogene 1; ADP-ribosylation factor-like protein 17		

For research use only.

GENE DESCRIPTION	ADP-ribosylation factor-like 17
GENE ID	641522
ORIGIN SPECIES	Human
ACCESSIONS	BC030570
VECTOR	pDONR223
NCBI SYMBOL	ARL17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggattcttattgagtttgatacagctggaaaaccaccatctgtataaattgaagctgggggagactgtgcctgccctacagtaggtttctgtgtggagacagtagaatataaaataaaccttcgctgtctgggatgtggcagccactcaaaatcagacctgtggcag
catTTTTCCAGAACACAAAAGGTGCCAGAAAGCCAGGAAGCACACATCAAGGCTCACTGCCAGCGGGGTGCTGCCAATAAATAGTACACGTGGAAATTTGGAATGTGGAAGGAGGTAGAAGTCATCTTCTCCCCATAGCAGCAGGTGTGAGGCTCTGTGTCAGCTGGACTC
catactccccaccagtcaccagcctggggaccgtggggctgcaaggacctcagcagcggtttcccaagtctcctgacttctccatcctctggaatcagctgtggtaaagtag