



Reference: PTXBC030601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PUS4; TruB pseudouridine (psi) synthase family member 1; TruB pseudouridine (psi) synthase homolog 1; TruB pseudouridine synthase family member 1		

For research use only.

GENE DESCRIPTION	TruB pseudouridine (psi) synthase homolog 1 (E. coli)
GENE ID	142940
ORIGIN SPECIES	Human
ACCESSIONS	BC030601
VECTOR	pDONR223
NCBI SYMBOL	TRUB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1