

ANKRD29-ankyrin repeat domain 29 Gene

Reference: PTXBC030622

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 29; ankyrin repeat domain 29		
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 29
GENE ID	147463
ORIGIN SPECIES	Human
ACCESSIONS	BC030622
VECTOR	pDONR223
NCBI SYMBOL	ANKRD29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcaggatgtcttcaagaaggaaactccactgccaatgctgcatctctgggcagccaggagaggaacctggcgctgctgaagctgctgttgaacagcggcgggtggacgtggactgcagagacagccatggcaccacactcctgatgtgtgctgctacgtggccacataga
ctgtgtgagggaactgtgtctgcaaggagcagacatcaatctccagagagagtcagggtacaactgcccattcttggcccccagcagccataatgatgtcgtgagatttctcttggatttggagcaccactgaatttaggacaaagacgagggcaccgcccgttggctgccagtc
gtacgggcacatgcaggtgtgtgagaccttctgaagcagggagcaaacatccatgaccaactttatgatggagccactgcccctcttctagctgccaaggtgtgttacttggatgttattcgattactgctggtctcagggagcaaaagtcaaccagccaaggcaggacgggacagcggccc
ctgtgtgatcgcgtccagatgggcccacagcaggtgtgtgctgctgctgcccggagccgacggcagcgtgctgctgacgacgatggcacaacagcattattgaaagcagccaacaaagggtataatgatgtcataaaagagttgcttaaatctcaccactcttggattttgaaga
atgggacatcagcgtccatgcagcagtgctcagtggaacattaaaacagttgctgctcctagaagcagggcgagaccatccctgagaacaaggccaatgaactccggcagaactaaccaaaaatgaactatattgctgctcctgagaagtaagaaggtcccagaagagc
taa