

## PRKD3-protein kinase D3 Gene

Reference: PTXBC030706

|                |                                                                                                                                          |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | EPK2; PKC-NU; PKD3; PRKCN; nPKC-NU; serine/threonine-protein kinase D3; protein kinase C nu type; protein kinase EPK2; protein kinase D3 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | protein kinase D3                                                                                |
| <b>GENE ID</b>          | 23683                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC030706                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | PRKD3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtctgcaataaattccctccatcagccagaagtctgtattaccacagctattctgtctgtctccagctcttccgtgttcaagtcctaaagacgggactctctgccgactctctaatggaagcttcagtcaccatcactcaccaactccagaggtcagtcacatcagtttcttctctctgcaaaattggcctcacacgggagagtggtaccattgaagccagggaactgtctttatctgtctgcaaggtctgtgtgtctccatggtttatcaaaaatttccagagtggtggtattcttggcatgtatgacaaaattcttctcttccatgacatgaactcagaaaacatttgcagctgattacctc  
aacagatgaatacatgaagagagacttagtggaagtggtcttctcagctttagccacagtagaagacttccagattcgtccacatactctctatgtacattcttacaagctctctctctgtgattactgtggtgagatgctctggggatttggtacgtcaaggactgaaatgtgaaggtgtggat  
taattaccataaacgatgtgcttcaagattccaaataactgtgtggaagtaagaaagagacgtctgtcaaatgtatctttaccaggacccggcctcagttccaagaccctacagccctgaatatgtagccctccagtggaagagtcacatgtccaccagggaaccaagtaagagaattcct  
tcttggagtggtcgcctaatctggaatggaaaagatgtaattgtcagagtgaaagttccacacacattgtctgtcactcttacacccgtccacgatgtcagtagtcaagcggttactgaaagccctcttcccaaggaatgcagtgtaagattgcaaatcaactgccataaacctgtg  
gcatcaaaagtaccaagagactgccttggaaggttacttcaatggagaaccttccagctctgggaacagatacagataccaatggatattgacaataatgacataaatagtgatagtagtgggtttggatgacacagaagagccatcacccccaagataagatgttctcttgatc  
catctgatctcgtatgtggaagagatgaagaagccgttaaaacaatcagtcacaaagcaataatattccgctaaggggtgtgtacaaatccatcaagcacacaaaggaagagcagcacaatggtggaaggaaggtgtggtgtccattacaccagcagggaataacctgagaaaga  
ggcattattggagacttgacagcaaatgtctaactatttcagaatgaatctggatcaaaagtattataaggaattccacttcagaattctccgcatacttccaccagagattcacaacatttcacaagcagcaatccacactgttttgaatcattactgatactatggtatactctgttggtg  
agaacaatggggacagctctcataatcctgttctgtgctgccactgagttgacttgatgtacacagagctgggaaaaagcaattcgcgaagccctcatgctgttactcctcaagcaagtgttgcattctccaggcgagggaagatcacaaagattgtctacaagtatctctgtatcta  
attgtcagattcaggagaatgtgatatcagttcttaccagatcttgcagatgaggtgtgtgttcaggccagtttggcatctgttatggagtaagcaactacagccgtttgcatattgtacacatttttaaaactggaaaatgtaa