

TMEM74-transmembrane protein 74 Gene

Reference: PTXBC030710

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET36; transmembrane protein 74		
For research use only.			

GENE DESCRIPTION	transmembrane protein 74
GENE ID	157753
ORIGIN SPECIES	Human
ACCESSIONS	BC030710
VECTOR	pDONR223
NCBI SYMBOL	TMEM74
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctccactacctgctaaagagcaaccaggcagacctctgtgatgccaggactggagttcaagaggctgcctggtgaccaggcagatacagcagccacaagagctgtctctgtgtcagaaacagtgatccacccaagagcaaccgagatggaagggtctaaac
ttagttcttccagcatccccctctctctctgcaaaacagtactctcagccagatgcttccaccaggacttccactcagggaacaaccaataacagcggaaaggaaagctgttaactgtcagccaggaaatagaaactctttacctatgtggacaaaacatcaactggagc
agcgaaccggagctgccatcagcaaaaggcataatcacctggggagcttggctgggaaatccaaatgagtggtcccaaggagctgccatcttggatctgaagaggagatgatacaagttcagaagccagctctcagggaagtctatagactatggtttcatcagcgcct
ctgttcttggctcactgggatcctgtcgtgatcatcttcaatcgtccacgggaagtgactgtgaccccaacactgtggcagccgggagatggagcgcctggagaaggagatgagggctggggctcactggaccgctgtgtgattgctggggctctgcctctcacgtggg
ggcgctcatcctgtcctgttgaatgatgccatgtggaagggggagctctatctgcaaacagatttcctcttccaaagagctgtcgaactctatggttcttcaactcaggatgaaaccagcacgaatgaaacactctggaactgtccttggtagaggagatgcgtgtgtaca
gagttaa