

CTD-2514C3.1-hypothetical LOC100134868 Gene

Reference: PTXBC030737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC100134868		

For research use only.

GENE DESCRIPTION	hypothetical LOC100134868
GENE ID	100134868
ORIGIN SPECIES	Human
ACCESSIONS	BC030737
VECTOR	pDONR223
NCBI SYMBOL	CTD-2514C3.1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagggttggcatgttccaggatctgtctatcgacttctcaggaggaatggaggcgctggacactgtcagaaggacttatag