

LOC133874-hypothetical gene LOC133874 Gene

Reference: PTXBC030767

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 5 open reading frame 58		

For research use only.

GENE DESCRIPTION	hypothetical gene LOC133874
GENE ID	133874
ORIGIN SPECIES	Human
ACCESSIONS	BC030767
VECTOR	pDONR223
NCBI SYMBOL	LOC133874
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagaaattgcagaaattccctgctcaggaaaaggtagaggcaaggattgacttaagatgggtaagaagcgtgttactgatcataagctaaatgtggacaaagtaattaaaaattacacaatttctcggagttgaagaagataaaagttaattga