



Reference: PTXBC030777

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UBP-M; UBPM; ubiquitin carboxyl-terminal hydrolase 16; deubiquitinating enzyme 16; ubiquitin specific protease 16; ubiquitin thioesterase 16; ubiquitin thiolesterase 16; ubiquitin-processing protease UB-P; ubiquitin-specific processing protease 16; ubiq		

For research use only.

GENE DESCRIPTION	ubiquitin specific peptidase 16
GENE ID	10600
ORIGIN SPECIES	Human
ACCESSIONS	BC030777
VECTOR	pDONR223
NCBI SYMBOL	USP16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]