



Reference: PTXBC030792

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDK5P35; CDK5R; NCK5A; p23; p25; p35; p35nck5a; cyclin-dependent kinase 5 activator 1; CDK5 activator 1; TPKII regulatory subunit; neuronal CDK5 activator; regulatory partner for CDK5 kinase; tau protein kinase II 23kDa subunit; cyclin dependent kinase 5		

GENE DESCRIPTION	cyclin-dependent kinase 5, regulatory subunit 1 (p35)
GENE ID	8851
ORIGIN SPECIES	Human
ACCESSIONS	BC030792
VECTOR	pDONR223
NCBI SYMBOL	CDK5R1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]