



Reference: PTXBC030825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf56; FTS and Hook-interacting protein; FHIP; family with sequence similarity 160 member A2		

For research use only.

GENE DESCRIPTION	family with sequence similarity 160, member A2
GENE ID	84067
ORIGIN SPECIES	Human
ACCESSIONS	BC030825
VECTOR	pDONR223
NCBI SYMBOL	FAM160A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]