

## SAP18-Sin3A-associated protein, 18kDa Gene

Reference: PTXBC030836

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | histone deacetylase complex subunit SAP18; histone deacetylase complex subunit SAP18; 2HOR0202; SAP18P; 18 kDa Sin3-associated polypeptide; Sin3A-associated protein, 18kDa; cell growth inhibiting protein 38; cell growth-inhibiting gene 38 protein; sin3-as |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Sin3A-associated protein, 18kDa                                                                   |
| <b>GENE ID</b>          | 10284                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC030836                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SAP18                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcggtagtgacgcgttaccagaggagaaattaagaaggagccagagaaccgacgcgcgagaagacatgccactgtgtacgggtcttcaccaccaataacggccgccaccacgaatggacgagttctccggggaaatgtaccgtccagcagtgatgcagatctac  
acttgatggatgcaacttgaagaactgacaagcttagtaaaagaagtctaccagaagctagaagaaggcactcactcaattttgcaatggttttacagatgttaaagacgtggctatcgagttaaggaaftggcagcaccatgtctggcagaaggggactgatgccatgac  
cctgcagtcgagaagttccagataggagattacttgacatagcaattaccctccaatcgggaccacactctcaggcgcatgagaccatattaa