

GPN3-GPN-loop GTPase 3 Gene

Reference: PTXBC031024

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPBD1C; GPN-loop GTPase 3; ATP-binding domain 1 family member C; protein x 0004		

For research use only.

GENE DESCRIPTION	GPN-loop GTPase 3
GENE ID	51184
ORIGIN SPECIES	Human
ACCESSIONS	BC031024
VECTOR	pDONR223
NCBI SYMBOL	GPN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcgggtatgcgcagctgggtcggcccgccggcagcgggaagagcacctactgtgccacatggtccagcactgtgaagccctcaaccgggtctgccaagtgttaaacctggatccagcagcagaacactcaactactcgtgatggtgacatccgggaactgacgagtg
gatgatgaatggagatgattctctgcgattcgggtcccaacggaggatgggtatgttcatggagtagtcttgcgaataatttgactggctggagaactgtctggccatgtaggagacgactatattctttgattggtccaggtcagattgattgtacactcacctgcctgtgatgaacagctg
gtccagcagctcgcagcagtgaggatccgagctctgtggagttttctgttgcattcagttcatgggtggagtcattcaagtttattctggcatctggcagccctgagtgccatgatctctctagaaatccgcaagtcacatcatgacaaaaatggatctgctgagtaaaaaagcaaaaagga
aattgagaaattttagatccagacatgtattctttattagaagattctacaagtgacttaagaagcaaaaaatcaagaactgactaaagctatattgtggactgattgatgactacagcatggttcgattttaccttacgatcagtcagatgaagaagcatgaacattgattgcagcatattgatt
tgccattcaatatggagaagacctagaatttaagaaccaagggaacgtgaagatgagctcttctctatgtttgacgaatatttcaagaatgccaggatgaatga