

BTG4-B-cell translocation gene 4 Gene

Reference: PTXBC031045

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein BTG4; APRO3; PC3B; B-cell translocation gene 4; BTG family member 4; protein PC3b; BTG anti-proliferation factor 4		

For research use only.

GENE DESCRIPTION	B-cell translocation gene 4
GENE ID	54766
ORIGIN SPECIES	Human
ACCESSIONS	BC031045
VECTOR	pDONR223
NCBI SYMBOL	BTG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagatgaattgcacaacagttttctgtcacaagattggtgaaaaacatgataaactaagtaaacagcaaatagaagactttgcagaaaagctgatgacgatctgtttgaaacatacagaagtcactggcactctgattgcccttctaaaggcgaagccttcaggtgcattcaggata
aacaacaatcagaataaagatccattctagaaggcagtggtgaaagttaattgttagattttctcacttgggacttccgaaggagatgaccatagggtagatccctttgaagtatgctgtaggtatggtgagaaaaacatccatttacagttgctcttttaaaggcagatggggaatgg
gaactatatcaacaatcagttatgcagtttagagcctcatcagacgtttccttggcacttctcgatgaagaaagttagcaaggaaacctgtgtcattcctaaagtcagcaatccgaagagtattatcaggtcaagagtggtccagttctttctatctatcttaattctaaaaag
aatgcactcattatgaaaaaagaagcaaaaaacatggaagaacaaactgtag