



Reference: PTXBC031051

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDH14; CDH14L; CDH24; cadherin-18; cadherin 18, type 2; cadherin-14; cadherin-like 24; ey-cadherin; cadherin 18		

For research use only.

GENE DESCRIPTION	cadherin 18, type 2
GENE ID	1016
ORIGIN SPECIES	Human
ACCESSIONS	BC031051
VECTOR	pDONR223
NCBI SYMBOL	CDH18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]