

LOC554207-hypothetical LOC554207 Gene

Reference: PTXBC031469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC554207
GENE ID	554207
ORIGIN SPECIES	Human
ACCESSIONS	BC031469
VECTOR	pDONR223
NCBI SYMBOL	LOC554207
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcaccagattctgttggctcaagaagcattcaagcttcaccccttcgtatttatgaatccggaggtactggagcagcgtgtatgtccagcatctggcattgttcaatcaagatgttag