



Reference: PTXBC031550

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein LMBR1L; LIMR; limb region 1 homolog-like; limb region 1 protein homolog-like; limb region 1-like protein -like; lipocalin-1 interacting membrane receptor; lipocalin-interacting membrane receptor; limb development membrane protein 1 like		

For research use only.

GENE DESCRIPTION	limb region 1 homolog (mouse)-like
GENE ID	55716
ORIGIN SPECIES	Human
ACCESSIONS	BC031550
VECTOR	pDONR223
NCBI SYMBOL	LMBR1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]