

C1orf183-chromosome 1 open reading frame 183 Gene

Reference: PTXBC031558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf183; protein FAM212B; family with sequence similarity 212 member B		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 183
GENE ID	55924
ORIGIN SPECIES	Human
ACCESSIONS	BC031558
VECTOR	pDONR223
NCBI SYMBOL	C1orf183
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgatggagcagggaatggactctatctccgtcgctcaaacaggagctgatgtccatgaaggaggtgggtgatggcttacaggatcagatgaactgcacatgaggtgcactgcaagaactgaagctcctcaggtgcagacagcactggaacagctggagatctctggagg
gggtctctgtccaggcagccctgaaggtccaggaccagtgccgagcaccctgtgtggagggtggcagaggtctccaggcccaagctctgttccccctccagtcacacttctctggcagcagcaccagttccatcccataggaagtgtctgtggaaggatagccccctgtccc
aggacacagccacatcaaaagctgtctcagcaggggccaagagcagtggaacggatgactggacctccacgttgatgtccggggccggaatcgacagcctctgtgttaggggacaacgttttgcagacctgtgtggcaattggctagacttgcagaactggagaagggtggg
gagaagggtgagactggggggcagctgaaccgaaggagagaaaggccagcccaaggagctggccgcaggtttgccctgacagcaaacatctttaagaagtcttgcgtagtgtgcggcctgaccgtgaccggctgtaaggagaagccaggctgggtgacacccatgggtcc
ctgagtcgccgaaccggcgtcacagaaggtaagaagcgagcctttccaaggctctgtgtcatttcccttcccaggccaccggggagcagagcgaggggagaatccccccacaagctgcccaaggccctggagcactcacctcaggattgatattaacacagctgtttgggtc
tga