

FAM27L-family with sequence similarity 27-like Gene

Reference: PTXBC031617

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	family with sequence similarity 27-like		

For research use only.

GENE DESCRIPTION	family with sequence similarity 27-like
GENE ID	284123
ORIGIN SPECIES	Human
ACCESSIONS	BC031617
VECTOR	pDONR223
NCBI SYMBOL	FAM27L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgttcgggatcaagatgaccacactccagccaaggataaaagcccccacaggagctcactgtcctgcaggagaggagcaggccacgtccaagaagatggctgtatgttttcacggctcttctctgagaagtgaagccacaccacagatacagcttgaagagggaagccgggaatgg
gagatggcaacaagccctgtcatgtgtgcctctctga