

GAGE12I-G antigen 12I Gene

Reference: PTXBC031628

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AL4; CT4.7; GAGE-7B; GAGE7B; G antigen 12I; G antigen 7B; cancer/testis antigen 4.7; g antigen 12F/G/I		

For research use only.

GENE DESCRIPTION	G antigen 12I
GENE ID	26748
ORIGIN SPECIES	Human
ACCESSIONS	BC031628
VECTOR	pDONR223
NCBI SYMBOL	GAGE12I
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgaggagatcgacattattggcctagaccaggcctgtacagcctctgaaatgattggcctatgcccagcaggtcagtgatgaagtggaccagcaacacctgaagaaggggaaccagcaactcaactcaggtcctgcagctgctcagggagaggatgagggagcctcaggtcaagggccgaagcctgaagctcatagcccgaacagggtcacccacagactgggtgtgagtgtaagatggtcctgatggcaggagatggaccgccaaatccagagggtgaaaacgcctgaagaaggtgaaaagcaatcacagtgtta