

SLC25A41-solute carrier family 25, member 41 Gene

Reference: PTXBC031671

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APC4; solute carrier family 25 member 41		
For research use only.			

GENE DESCRIPTION	solute carrier family 25, member 41
GENE ID	284427
ORIGIN SPECIES	Human
ACCESSIONS	BC031671
VECTOR	pDONR223
NCBI SYMBOL	SLC25A41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggctctactctccaagacgaacttcaccaactgctggggggctacagagcatgtccaggaggcggtctccctgtggcggggcaacggcatcaacgtgctcaagattgctctgagtagccatcaagttctccgtattcgagcagtgcaagaattactctgtggaat
acaagggtcccccgccttcaggagcgtctctgtggtccctgtgtggccatctcccagacctcatcaacccatggagggtgctgaagacggtgtgaccttgcgtcggacgggcccagtaagggtgctgactgcgccaggcagatctgcagcgagaggcaccgc
gccccttaccggtgctacgtcccaatatgtcggcatcatccctatgctgcacgacctggctgtatgagatgctccaggtctctgtggtgaagtgcagggagatgggggacccagtgccctggctcagctgtgtgtgacgctatccagacctgtggccagatggccag
ctacccactgactctgtgtgcaccaggatgcaggccaagataccgtggagggtcaaatcccacatgcgggagtcctccagcgatcctggccagcagggtggctagggtgtaccgagggatgacccccacgctactgaagggtttaccagcaggtggcatcagctatgtg
gtgtacgaagccatgaagaaacctgggcatatag