

LOC643406-hypothetical protein LOC643406 Gene

Reference: PTXBC031676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC643406		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC643406
GENE ID	643406
ORIGIN SPECIES	Human
ACCESSIONS	BC031676
VECTOR	pDONR223
NCBI SYMBOL	LOC643406
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgcctttgcactgtccatgctctcagccaccacccctccccctgccattgcagtggctaacactggagacgaagaccaagacaccgcctttctccagctacctctcaaatcagcacaggcaaggacaaggcctcaatccacaattgtgaagatggaccctggcccatgg
gatggtcagacacgcctgccagctatctgcaggcgaagaggtcagaagaggttaggggcctgaaggacatctgtctccatgtccatagcaggctatttctgtccatga