

RFPL3-ret finger protein-like 3 Gene

Reference: PTXBC031689

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ret finger protein-like 3; ret finger protein like 3		
For research use only.			

GENE DESCRIPTION	ret finger protein-like 3
GENE ID	10738
ORIGIN SPECIES	Human
ACCESSIONS	BC031689
VECTOR	pDONR223
NCBI SYMBOL	RFPL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgactcttccaagaagcaagcagctgtcccgctgtcagactatctggaaaaaccaatgtccctggagtgatgcacgtctgctcaagtgcacaaatcgtgcagaaggagcccatggggaggatctgtttgtgtgtgttccatggtctcagaggaaacaaatca
ggcccaatcgccagctagagagctgtgttccacatcaaggaactggagcccaagatgaagaagattctacagatgaaccaagatcggaagtccaagtggatagacctggatgccacacagccaaacttctctcatttctgacacctcaggagcgtccgaagtgggc
tcacacacagaatcggaagacctgtccgagagatttgacgtgtccgtttgcatctgggtccctcgtttacctgtggccgccactactggaggtggacgtgggaacaagcacagaatgggacctggagctgtgcagagaatctgttactgcaagggaagatccagctgaccac
agagcttggattctggactgtgagtttgaggagtggaagccgctctcgtccagcaggtgccgctgacttctcttagtagaccgcaagttacagcagtgagggttttttggatattggcatggcagcagaacgtttctttttgatgctgaaagtgttccatgtctatactcaggagcgtct
ctgctgaggagccactgcgccattttggctcttcaattccacctaattggtgatcagggtgtcttgagcatctgtcctttgatgaactcaggcactactgatctccagtcctgctgggagggccaaataa