



Reference: PTXBC031828

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF2A; CDCBM3; HK2; KIF2; Kinesin, heavy chain, 2; kinesin heavy chain member 2A; kinesin-2; kinesin family member 2A		

For research use only.

GENE DESCRIPTION	kinesin heavy chain member 2A
GENE ID	3796
ORIGIN SPECIES	Human
ACCESSIONS	BC031828
VECTOR	pDONR223
NCBI SYMBOL	KIF2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]