

C15orf50-chromosome 15 open reading frame 50 Gene

Reference: PTXBC031958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf50; long intergenic non-protein coding RNA 593		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 50
GENE ID	414926
ORIGIN SPECIES	Human
ACCESSIONS	BC031958
VECTOR	pDONR223
NCBI SYMBOL	C15orf50
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgataagctctgggtggagacgcctatctcagaggagctccaaacctgctcaaggcctttctgcccctctgcaactttctgcccaagcaattagcctcagggagaagaggaaattggaacacggacacgcacagaggaaagaccaagtga