



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MRCL3-myosin regulatory light chain MRCL3 Gene

Reference: PTXBC031972

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRCL3; HEL-S-24; MLC-2B; MLCB; MRLC3; MYL2B; myosin regulatory light chain 12A; epididymis secretory protein Li 24; myosin RLC; myosin regulatory light chain 2, nonsarcomeric; myosin regulatory light chain 3; myosin regulatory light chain MRLC3; myosin, 1		

For research use only.

GENE DESCRIPTION	myosin regulatory light chain MRCL3
GENE ID	10627
ORIGIN SPECIES	Human
ACCESSIONS	BC031972
VECTOR	pDONR223
NCBI SYMBOL	MRCL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcagcaaaaagcaaaagaccaaagaccaaagaagcgccctcagcgtgcacaatccaatgtgttgcctatgtttgaccatgcgcagatcaggaggttcaaaagaggccctcaacatgattgatcaaacagagatggtttcatcgcagaagaaagatttcgatgatgcttgcctcattggggaagaattccaagttagatgattgatgccattgaatgagatgaggcgcgccatcaattccaatgtgttccatcagatgttggggaaagttaaaggcgacagctctggaatgatgtcatcagaagaatgcttgcctgttccttgatgaagaagcaactggtcccatcatcaggaagattactctagagacgtctcacaacacatctgctcgtctttacagatgcagaaggtggatgagctctacagagaagcaccattatgataaaaagggttgaatcattcaatcagtgctcagctctgaaccatctgaaccatctgaccaaagacaagaatgactga

More info: Please contact us for any question: contact@proteogenix.fr