

THRSP-thyroid hormone responsive (SPOT14 homolog, rat) Gene

Reference: PTXBC031989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LPGP1; Lpgp; S14; SPOT14; THRP; thyroid hormone-inducible hepatic protein; SPOT14 homolog; lipogenic protein 1; spot 14 protein; thyroid hormone responsive (SPOT14 homolog, rat); thyroid hormone responsive SPOT14; thyroid hormone responsive		

For research use only.

GENE DESCRIPTION	thyroid hormone responsive (SPOT14 homolog, rat)
GENE ID	7069
ORIGIN SPECIES	Human
ACCESSIONS	BC031989
VECTOR	pDONR223
NCBI SYMBOL	THRSP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagggtgtaaccaagcgttacccaagaactgctgctgacgctcatggaccggtatgcagccgaggtgcacaacatggagcaggtggtgatccccagcctctgcgggacgtgcagctgagtgaggcctggggccaggccaggtgagggccctgatctctacactactt
caccatgctcaaggccatctgtgtggaatgtagacatgggctgctgcccggggaggagtgagcagccaaggtggcaggcagcgaagagaatggaaaccgagagacagaggagtcgaggacgagagtgccctcaggagagctggacctggaagccagttccactgcacttctcc
agcctccatcacatctcatgcacctcaccgagaagccaggaggtgacaaggaaataccaggaaatgacgggacaagtttgtag