

C22orf23-chromosome 22 open reading frame 23 Gene

Reference: PTXBC031998

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EVG1; dJ1039K5.6; UPF0193 protein EVG1; chromosome 22 open reading frame 23		
For research use only.			

GENE DESCRIPTION	chromosome 22 open reading frame 23
GENE ID	84645
ORIGIN SPECIES	Human
ACCESSIONS	BC031998
VECTOR	pDONR223
NCBI SYMBOL	C22orf23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcacagaagcagatggagtagtgaccaagggaactgggtccggcgcccccagaccatctacacccggggacctgcgagctgctcagatgatgatgaaggaaatccaaactgacgaacatccagcagcgcacatcatggacatcatgaaaagaggagatgctt
tgccccctacagtgcagcccaatccagccagagagtgcttacctccaagcaaatagcctgcccctctacctgcccctctctcgcagcccgctccacctccggcctgccaacatgtgtcaagccaatgggcctacagccgggagcagttcaagcctcaagccaccagggaattgg
agaaggagaaacaaagactccaaatatcttggccacggggaaggacatggaggaacgaaaaaaggccccctctgcacgacagaaggctccagccctgagctagaccgatttgaagagctggtgaaggaaatccagagaggaaagaattctgctgacatggagccctg
ggacagggcaaacagtagggaatcatcttctgaaatctccagaaactccgggaatggaagacattgaccacagaaggatgagggaacttaggaagggtcttggccaccactaa