

FAM27E3-family with sequence similarity 27, member E3 Gene

Reference: PTXBC032035

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	family with sequence similarity 27, member E3		

For research use only.

GENE DESCRIPTION	family with sequence similarity 27, member E3
GENE ID	100131997
ORIGIN SPECIES	Human
ACCESSIONS	BC032035
VECTOR	pDONR223
NCBI SYMBOL	FAM27E3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggatattcaattgctcgggaccgaagaatctctcactgtggccaggccttcacacacccaaagcggaaccgctcggcgaaaaggattgacaaccgtctcatgacccaggcagagaggcagaacaggctcaccaagacaggccgcatgcgagaaacagctttgtgg
tgcacagggcacattcggccaaggacacacacgcacacgggcacacacacacaaaccgacagagaggggaagaacacacagagactgagagacagagagagaagagagaatgggagacacacacacatacacacacagacacacacagagctttatag