

ZNF44-zinc finger protein 44 Gene

Reference: PTXBC032246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GIOT-2; KOX7; ZNF; ZNF504; ZNF55; ZNF58; zinc finger protein 44; gonadotropin inducible transcription repressor-2; gonadotropin-inducible ovary transcription repressor 2; zinc finger protein 55; zinc finger protein 58; zinc finger protein KOX7; zinc finge		

For research use only.

GENE DESCRIPTION	zinc finger protein 44
GENE ID	51710
ORIGIN SPECIES	Human
ACCESSIONS	BC032246
VECTOR	pDONR223
NCBI SYMBOL	ZNF44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtgttatggaactttctgggctaccctaagatgctggaagctgccaatctcatggaggccctagtgatcgcgccttgggtcactctccagaggacagcctgaggtgttgagtgaggcctccaaaggatcaggactcagtgcccttgaggatgtgctgtgaactc
acccatgaggagtgaggctttgctgggtccatcacagaagaatctctacagagatgtgatgcgagaaccattaggaaacctgaactgtataggaatgaaatgggaaaccagaacattgatgacgacacaaatctcaggagaaatccaagggttatcgagacatggttacagaattcatt
ttgatccattatggtcacctgtagctttgtctgatggaggtttgttacagttttcaacaggacaggtctccctgtcaccagcgtggagttcagtag