

BATF-basic leucine zipper transcription factor, ATF-like Gene

Reference: PTXBC032294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B-ATF; BATF1; SFA-2; SFA2; basic leucine zipper transcriptional factor ATF-like; B-cell-activating transcription factor; SF-HT-activated gene 2 protein; activating transcription factor B; basic leucine zipper transcription factor, ATF-like; basic leucine		

For research use only.

GENE DESCRIPTION	basic leucine zipper transcription factor, ATF-like
GENE ID	10538
ORIGIN SPECIES	Human
ACCESSIONS	BC032294
VECTOR	pDONR223
NCBI SYMBOL	BATF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctcacagctccgacagcagtgactccagctcagccgctctctccccctggcaaacaggactcatctgatgatgtgagaagagttcagaggaggagaaaaatctattgccgcccaagaagagccgacagaggcagacacagaaggccgacacacctgcacctggagagcgaa
gacctggagaacacagaacgcggctctacgcaaggagatcaagcagctcacagagggaactgaagtacttcacgtcgggtgctgaacagccacgagccctgtgctcgggtgctgcccagcacgcccctgccccccgaggtggtgtacagcgccacgcattccaccaacctcatgtc
agctccccgcgcttcagccctga