

CISD2-CDGSH iron sulfur domain 2 Gene

Reference: PTXBC032300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERIS; Miner1; NAF-1; WFS2; ZCD2; CDGSH iron-sulfur domain-containing protein 2; endoplasmic reticulum intermembrane small protein; mitoNEET-related 1 protein; nutrient-deprivation autophagy factor-1; zinc finger, CDGSH-type domain 2; CDGSH iron sulfur dom		

For research use only.

GENE DESCRIPTION	CDGSH iron sulfur domain 2
GENE ID	493856
ORIGIN SPECIES	Human
ACCESSIONS	BC032300
VECTOR	pDONR223
NCBI SYMBOL	CISD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgctggagagcgtggcccgatcgtgaaggtgcagtcctgcatactgaagcggcctccagtcctgaaagcattaccgggttcgtaggtcacagttcagaatggcttcggttattgccttccttggtgtactcgcactcttggtctacctgcagttcgtccattcctccgaaga
agaacaacagaaggatagcttgattaatcttaaaatacaaaaggaaaatccgaaagttagtgaatgaataaacattgaagatttgtcttactaaagcagcttatttaggtgttggcgttctaaaacgtttcctgcctgcgatgttcacataataacacaatgaattgacaggagataatg
gggtccactaatactgaagaagaagaagtataa