

SPIC-Spi-C transcription factor (Spi-1/PU.1 related) Gene

Reference: PTXBC032317

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPI-C; transcription factor Spi-C; Spi-C transcription factor (Spi-1/PU.1 related); Spi-C transcription factor		

For research use only.

GENE DESCRIPTION	Spi-C transcription factor (Spi-1/PU.1 related)
GENE ID	121599
ORIGIN SPECIES	Human
ACCESSIONS	BC032317
VECTOR	pDONR223
NCBI SYMBOL	SPIC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgtgtgtgaacaagaagctgggtcaagcattgaagatgctttgaggtctgaggcaacattcaactggagatcttcagtactcgccagattacagaattacctggctttaatcaaccatgctcctcatgtcaaggaaattccagctgctatggagtgtgcctacagaggagcctgtctataattggagaacggttaataacagtgctgctgactctattttgaaggaaatattcatcaatctctgcagaacataactgaaaaccagctgggtacacccactctctccagcaaaaggggggaaagcgaggagaagctccgactgtttgaataccttcacgaatccctgtataatccggagatggcatctgtattcagtggttagataaaaccaaaggcatcttcagttgtatcaaaaaacaagaaaaacttgcggagctttgggggaaaagaaaggcaacaggaaagacatgacttaccagaaaatggccagggcactcagaattacggaagaagtggggaaattacaaatccggaggaaagctgacttaccagttcagtgaggccattctcaagactctccatcctatttctggggaagagatcttctattcagtggtttcaacctgatcaagaatatctcagtttaataaactggaatgcaaatataattacatatgccaatccatgagctaaatccatgatgtctaa