

C12orf45-chromosome 12 open reading frame 45 Gene

Reference: PTXBC032326

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C12orf45; chromosome 12 open reading frame 45		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 45
GENE ID	121053
ORIGIN SPECIES	Human
ACCESSIONS	BC032326
VECTOR	pDONR223
NCBI SYMBOL	C12orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtccatggcaagccccaggttagcccgagttgttcgtcgcacccgggattctcaggagtcctcaggtgtccaaggagctgctgacggcggaagcgcggcgagggtatgggacaggttgcctcaactcccaactaagtcagaaagacctccactctcaaa
cagttcggatagagaggagtccttattggaccaggtacagacatttctccacagatggcagggcgaatgaaagctaaagaaagaaatggcagctgcaccacctggtcgtttcaatattgaaaacattgatggcctcatagtaaaagtatacaaatggatgtgctttgttgatgaat
cagtcggattcaaaaaggtggacaggttcagaagagaggttcacaagacagttcagagaaacagttcagaatcagaagacgaagatgacagcatcccatctgaagtcacatagataacattagctcccaattctgaaggtggaaaaggcaagattgaagttttggacagtcagcaagta
aaaaaaaaagaatagccaataa