



Reference: PTXBC032374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MCM6 minichromosome maintenance deficient 6 (MIS5 homolog, <i>S. pombe</i>); DNA replication licensing factor MCM6; MCG40308; Mis5; P105MCM; minichromosome maintenance deficient (<i>mis5</i> , <i>S. pombe</i>) 6; minichromosome maintenance complex component 6		

For research use only.

GENE DESCRIPTION	minichromosome maintenance complex component 6
GENE ID	4175
ORIGIN SPECIES	Human
ACCESSIONS	BC032374
VECTOR	pDONR223
NCBI SYMBOL	MCM6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]