

PLXNA1-plexin A1 Gene

Reference: PTXBC032432

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOV; NOVP; PLEXIN-A1; PLXN1; plexin 1; semaphorin receptor NOV; plexin A1		
For research use only.			

GENE DESCRIPTION	plexin A1
GENE ID	5361
ORIGIN SPECIES	Human
ACCESSIONS	BC032432
VECTOR	pDONR223
NCBI SYMBOL	PLXNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggggtgcaccacaggcctcatgaagcagttccacatgggcgtgtggcgtggccaccacagagcacatggcgtgtgtctagggcgcaagcacttttagcagtatctgtttacatgcgcaaggatcaagccgactacctgtgtgtctactgggacagcagttccgagcta
ctccgtacctccctctgccaggctgtggagttaggccccagtcctacttgcactgtgtccactgtgtctcctaactgtgcagcacctgggagctctggcctggggtggagccctggtag